Abb. I, Fig. I
Resin I

Abb. II, Fig. II
Resin II



Abb. III, Fig. III
Resin III

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(EN)

INSTRUCTIONS FOR USE for the OptiMed CibuFlex, Ösophagus Stent

Please Note! The product is sterile packaged (EO sterilisation), designed for single use and must not be re-sterilised. The product must be stored in a cool, dry place and protected from light. This product must not be used if either the product or packaging are damaged. **Only a doctor familiar with its application may use this product. This medical product is designed for single use only. It must not be reprocessed, cleaned, disinfected or sterilised!** After reprocessing, the product is no longer suitable for use and may constitute a biological risk for both patient and user: owing to the surface composition and design of the product, full effectiveness of cleaning after use cannot be guaranteed and contamination cannot be ruled out if the product is re-used. Re-sterilisation may impair the surface composition and material properties of the plastics used. In this case the original function of the product cannot be guaranteed.

The CibuFlex, Ösophagus Stent is a self-expanding Nitinol stent with fixed ends between which is an ePTFE-coated section. The system comprises two elements:
1. **CibuFlex, Ösophagus Stent**
2. **Application device 22F**

The physician performing the stenting procedure must be consulted if the patient is known to have a nickel allergy. Indications, contraindications and application of the above product will always comply with current medical standards. The guidelines and recommendations of the relevant specialist medical associations must be observed.

- Indications**
- Non-vascular indication
 - Dysphagia in malignant incurable oesophageal stenosis
 - Malignant oesophago-respiratory fistula

- Absolute Contraindications**
- Occlusions which cannot be rechannelled (guide wire cannot be advanced to /through stenosis)
 - Obstructions with proximal end localised less than 3 cm below the upper oesophageal sphincter
 - Application at oesophago-gastric junction

Relative Contraindications

- Application at the junction oesophago-gastric

Possible Complications
Only doctors familiar with the possible complications may apply this product. Complications may occur at any time during or following the procedure. Possible complications include:

- Impossibility of widening stenosis
- Perforation of the stent/oesophagus
- Acute occlusion of the oesophagus
- Oesophageal haemorrhage
- Re-dysphagia from tumour growth over the uncovered part of the stent / upper and lower edge
- Stent migration
- Incorrect positioning of the stent
- Inadequate anchoring /dislocation resulting from too small a stent diameter
- Aspiration
- Blood pressure dysregulation
- Cardiac arrhythmia
- Allergic reaction to contrast medium
- Infection/sepsis
- Death

- Please Note!**
- The lumen diameter and length of lesion must be measured prior to every stent implantation.
 - The stent length selected should generally be 4 cm (2.0 cm distal, 2.0 cm proximal) longer than the stenosis /lesion.
 - Depending on the morphology of the lumen, prior to stent implantation balloon dilation should be performed to facilitate passage of the application device.

- If balloon dilation is necessary after positioning the stent (i.e. waiting in cases of major stenosis) this is generally performed using X-ray monitoring until the optimum lumen width is reached or at maximum to the defined diameter of the stent in question.
- **Note:** Balloon dilation must be performed carefully and slowly (15-30 minutes) until the patient registers pain; rapid and complete balloon dilation is contraindicated and increases risk of perforation.

- The CibuFlex, Ösophagus Stent should not be implanted at points where other stents of other materials may touch the stent or are in its vicinity, as this may cause electrolytic corrosion.
- The CibuFlex, Ösophagus Stent is not designed for use in the central circulatory system (as per EC Directive 93/42/EC, Annex IX, Para 1.7.).
- During manipulation of the application device a guide wire must always be present in the guide wire lumen.
- Chemotherapy and radiotherapy carry the risk of stent migration as a result of tumour shrinkage, stent erosion and/or muscle haemorrhage.

Please refer to the OptiMed brochures or product labels for the precise details of the above products.
A physician and patient form is included in the stent packaging materials. Fill in the stent data, indication, size and area of the implantation. For this the removable label on the package can be used.

- Preparation**
- Check packaging for damage. If the packaging is damaged, the product must not be used.
 - Remove product from the package using aseptic methods and continue to handle aseptically.
 - Check product for damage. A damaged product must no longer be used.
 - Ensure that products not belonging to the set are compatible and read their instructions for use carefully.

Note
Do not link the application device!
When using the releasing stent, ensure that the flexible application device is positioned as straight as possible along its length. To straighten the catheter, slightly withdraw the proximal handle.

- Keep the handle of the application device steady while releasing the stent.
- Do not rotate or twist the application device while releasing the stent.
- A partially or fully expanded stent cannot be retrieved into the application device!

Handling
Preparation and introduction of the application device (Fig. 1)
1. Remove the closed application system from transport packaging, holding the Y-adapter (2) to do so.
2. Remove the transport carrier.
3. Flush the central lumen and Y-adapter with physiological saline solution (non-heparinised).
4. Introduce the application device, using the guide wire (see label for diameter of guide wire).
5. Monitor the positioning of the application device by X-ray and ensure that the platinum markers (5) overlap the obstruction proximally and distally.

Releasing the stent (Fig. II)
6. Fully unscrew the fastening screw (2a) by turning it several times counter-clockwise (it is not possible to rotate the screw through 360° as necessary).
7. Leave the lumen centreline and the application device in Y with a solution of saline physiological (non-heparinised).
8. Insert the system of application to long of the fit-guide (4) in a distal direction (5) in a distal direction towards the tip (6), monitoring the procedure by X-ray.
9. Further procedure according to current medical standards.

A partially released stent cannot be retrieved into the application device.

Removing the application device (Fig. III)
9. When the stent is fully released, close the application device. Move the sheath (3) in a distal direction until the Y-adapter (2) has reached its original position.
10. Now withdraw the application device via the guide wire in situ.
11. Perform an angiogram to check the correct positioning of the stent.

Repositioning the stent
If the stent has expanded too far below and too far away from the stricture after release, the extraction loop (4a) can be used to correct its position. To do this, use the operating channel of an endoscope to introduce a forceps/snare, grasp the extraction loop at the proximal end of the stent and carefully draw the stent into the desired position.

Withdrawing the stent
To withdraw the stent (4), use the operating channel of an endoscope to introduce a forceps / snare, grasp the extraction loop at the proximal end of the stent and carefully withdraw the endoscope with the stent and forceps / snare from the oesophagus.

Note: In the event of a complete, a repairable or potentially repairable incident, all products connected with the incident and made by OptiMed or any other brand must be preserved wherever possible. If required, they must be handed over to the responsible authorities or, if necessary, to the manufacturer itself for further analysis.

OptiMed products are not entitled to change the above terms and conditions, to extend the liability, or to assume additional obligations concerning the product.

OptiMed guarantees that its products are manufactured with the greatest possible care. THIS IS THE ONLY VALID GUARANTEE AND SUBSTITUTES ALL OTHER GUARANTEE DECLARATIONS WHICH HAVE BEEN GIVEN.
It should be noted that no product is always effective under all circumstances; this is due to the biological differences of the persons to be treated.

Components of OptiMed sets as well as individual OptiMed products are interchangeable, provided that information on sizes is observed. Before using individual OptiMed products or sets with non-OptiMed products, the user must ensure the compatibility of the individual products with the specific application.
OptiMed has no influence on how the product is used, the diagnosis of the patient, nor how the product is handled outside of the company. OptiMed cannot guarantee the effectiveness nor the absence of complications. Therefore, OptiMed assumes no liability for damages and costs. OptiMed will replace products which have a defect for which OptiMed is responsible. OptiMed cannot be held liable for consequential damages of any kind which was caused by reprocessing or reprocessing the product.

In the event of a complete, a repairable or potentially repairable incident, all products connected with the incident and made by OptiMed or any other brand must be preserved wherever possible. If required, they must be handed over to the responsible authorities or, if necessary, to the manufacturer itself for further analysis.

OptiMed products are not entitled to change the above terms and conditions, to extend the liability, or to assume additional obligations concerning the product.

Products are subject to alterations.

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