GERBACHSANWEISUNG

OptiMed® Ureterale, Uretrale Endoprothesen, Intra-Operative Ureterale, Ureterostomie- und Steurbare Ureterkatheter

Beachten Sie!

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- Ureterkathete
 - Distales Ende geschlossen, standard und Tiemann opening
 - Closed or open, via ureteroscope
 - Closed or open for pediatric use
 - Tamour stent, closed or open
 - MultiTechnique, MultiLength
 - Shearable ureteral stent „Dirigent Optimaic“ closed or open
 - For kidney transplants according to Pevring
 - Open
 - Closed or open and Tiemann
 - Closed or both ends
 - Ferri des davis and solon
- Ureterostomie (Endoprothesen)
 - Offen
 - Ureterostomie /offen und Tiemann
- Intra-operative Ureterkathete
 - Beidseitig geschlossen

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Intra-operative Ureterkathete

- Offene operative Eingriffe am Uter

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- Reine die Ureterale, die des Puhers und die Führungsdrähte mit physiologischer Kochsalzlösung.
- Identify the flexible end of the guide wire.
- Insert the flexible guide wire end through the working channel of the optical instrument used.
- Locate and probe the ureter and position the guide wire end in the renal pelvis until the distal pigtail is introduced in the ureteral stent through the postioned guide wire.
- Push the ureteral stent forward with the help of the guide wire, using the postioned guide wire, up to the renal pelvis until the distal pigtail is introduced in the ureteral stent.
- Retract the ureteral stent until the distal pigtail is in the renal pelvis.
- The proximal end of the ureteral stent (black band before last of the stent) must be placed in the renal pelvis.
- After position control providing the image (X-rays, ultrasound), remove guide wire and trigger the dropping mechanism, by pressing the application system base until it releases the ureteral stent.
- The cytoscope is subsequently removed.

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EN

INSTRUCTIONS FOR USE

OptiMed® Ureterale, Uretrale Endoprothesen, Intra-Operative Ureteral Stent, Ureterostomy Stent and Controllable Ureteral Stent

Please Note!

The product is sterile packaged (EO sterilization), designed for single use and must not be re-sterilized. The product must be stored in a cool, dry place and protected from light. This product must not be used if either product or packaging are damaged. Only a doctor familiar with urological instruments. The instruction may only be used by a doctor who is familiar with the product. This medical product is designed for single use only. It must not be reprocessed, cleaned, disinfected or sterilized! After reprocessing, the product is no longer suitable for use and may constitute a biological risk for both patient and user. Saving to the surface composition and design of the product, but effectiveness of cleaning after use cannot be guaranteed and contaminants cannot be ruled out if the product is re-used. Re-sterilization may impair the surface composition and material properties of the plastics used. In this case the original function of the product cannot be guaranteed.

Description/ Composition

- Ureteral Stent
 - Distal end closed, standard and Tiemann opening
 - Closed or open, via ureteroscope
 - Closed or open for pediatric use
 - Tamour stent, closed or open
 - MultiTechnique, MultiLength
 - Shearable ureteral stent „Dirigent Optimaic“ closed or open
 - For kidney transplants according to Pevring
 - Open
 - Closed or open and Tiemann
 - Closed or both ends
 - Ferri des davis and solon
- Ureterostomy (Endoprothesis)
 - Open
 - Ureterostomy /open and Tiemann
- Intra-operative Ureteral Stent
 - Beidseitig geschlossen

Note

If the pusher has radiopaque markings, it can be clearly identified under fluoroscopy and distinguished from the stent portion (then under an X-ray examination this serves towards the identification of what is being pushed and/or what is characteristic).

OptiMed Ureteral stents, uretrale endoprothesen, intra-operative ureteral stents and ureterostomy stents are made out of a highly developed optical instrument, which is known as X-ray. They have an extremely low tendency to encrustation, are biocompatible, have a high degree of uniaxial stability and a homogeneous surface.

The tumour stent is additionally strengthened with polyimide in the stent.

Biocompatibility tests of the OptiMed ureteral stents show a potential period of use of up to one year. However, the functional efficacy of the ureteral stents must be regularly checked in accordance with the guidelines and recommendations of specialist urological associations.

OptiMed Ureteral stents may remain in situ in the patient for a maximum of 29 days.

Indications, contraindications and application of the ureteral stent must always comply with current medical standards. The guidelines and recommendations of the relevant specialist medical associations must be observed.

Indications

- Ureteral Stent /Pediatric Ureteral Stent**
 - General obstructions and stenoses in the area of the ureter
 - Supportive measures in the therapy of renal calculi and ureteritis
 - After surgical interventions in the urogenital tract
 - Trauma

Tumour Ureteral Stent

- Palliative measure, for inoperable tumours with a malignant dignity. The same application procedure applies as for the conventional ureteral stent.

Urological Endoprothesis

- Indications like for the Ureterkathete, additionally however, as a measure for short-term drainage and diagnostic urine sampling

Intra-Operative Ureteral Stent

- Open operative interventions in the ureter

Ureterostomy Stent (Endoprothesis)

- Inzision of the Ureters, wobei die Anwendung wie bei der konventionellen Ureterkathete erfolgt. However, care must be taken to ensure that the sleeve is attached, on the same level with the incision

Ureteral Stent Set for Kidney Transplants According to Pevring

- Optimal percutaneous drainage after the kidney transplant

Contraindications

- No safe access possible
- Non-therapizable drainage reductions

Possible Complications

Only doctors familiar with the possible complications may apply this product. Complications may occur at any time during or following the procedure. Possible complications include:

- Encrustation /Obstruction of the ureteral stent
- Migration of the ureteral stent
- Intolerably caused infections /Septis
- Organ perforation /Obstruction
- Haematuria /Haematurie

Note: Products whose drainage effect is impaired must be changed. The statutory hygiene regulations must be adhered to.

Registration

Please refer to the OptiMed brochures or product labels for the precise details of the above products. The packaging consists of a doctor's and patient card where the exact details, indication and implantation areas are entered. To this end, the removable label on the package can also be used.

Phthalates

This medical product may contain phthalates. Please check the label on the packaging for information. Phthalates (plasticizers) are used to increase the flexibility of the material. The product has been designed and produced in such a way as to minimize risks arising from substances released by the product. This also applies to the use of the product with children and pregnant or breastfeeding women.

Preparation

- Measuring the ureter to determine the correct size and length of the product.
- Check packaging for damage if the packaging is damaged, the product must not be used.
- Remove product from the packaging using aseptic methods and continue to handle aseptically.
- Check product for damage. A damaged product must no longer be used.
- Ensure that products not belonging to the set are not compatible and read their instructions for use carefully.

Handling Ureteral Stent (Open Version)

- Reine die Ureterale, die des Puhers und die guide wire with physiological saline.
- Identify the flexible end of the guide wire.
- Insert the flexible guide wire end through the working channel of the optical instrument used.
- Locate and probe the ureter and position the guide wire end in the renal pelvis until the distal pigtail is introduced in the ureteral stent through the postioned guide wire.
- Push the ureteral stent forward with the help of the guide wire, using the postioned guide wire, up to the renal pelvis until the distal pigtail is introduced in the ureteral stent.
- Retract the ureteral stent until the distal pigtail is in the renal pelvis.
- The proximal end of the ureteral stent (black band before last of the stent) must be placed in the renal pelvis.
- After position control providing the image (X-rays, ultrasound), remove guide wire and trigger the dropping mechanism, by pressing the application system base until it releases the ureteral stent.
- The cytoscope is subsequently removed.

Handling Shearable Ureteral Stent „Dirigent Optimaic“ (Closed Version)

- Reine die U