


## Chiba Needle

The combination of a specially cut cannula and modified stylet guarantees excellent visibility of the tip under ultrasound.

available with bevel-edged Cut or US-Cut

Chiba Needle, US-Cut



Bevel-edged

Diameter mm	Diameter G	Length mm	Order-Nr.
0.7	22	150	1411-1150
0.7	22	220	1411-1220

Box: 10 PCS

US-Cut

Diameter mm	Diameter G	Length mm	Order-Nr.
0.7	22	90	1413-1090
0.7	22	150	1413-1150
0.7	22	220	1413-1220

Box: 10 PCS