

low friction  
easy entry

DLC

Diamond Like  
Carbon Coating

## CarboSoft

Ureteral Stent

Die Langzeit-Ureterschleife mit innovativer DLC-Beschichtung für maximale Patientenverträglichkeit durch verringerte Biofilmbildung

*The long-term ureteral stent with innovative DLC-coating for maximum patient tolerance due to reduced biofilm formation*

6CH

7CH


8CH

**OptiMed**  
global care



- DLC-Beschichtung (Diamond Like Carbon)
    - extrem widerstandsfähig
    - hohe Biokompatibilität
  - Ultraglatte Kohlenstoff beschichtete Oberfläche
    - reduzierter Reibungswiderstand
    - erlaubt ein schonendes Einführen der Ureterschiene
    - maximiert den Patientenkomfort und die Patientenverträglichkeit
  - Verringerung des Risikos von Infektionen, Biofilmbildung und Inkrustationen
  - Hervorragend geeignet für die Langzeitanwendung
- *DLC coating (Diamond Like Carbon)*
    - *extremely resistant*
    - *high biocompatibility*
  - *Ultrasmooth carbon coated surface*
    - *significantly reduced friction*
    - *enables smooth introduction*
    - *maximized patient comfort and tolerance*
  - *Reduced crystalline biofilm formation, risk of encrustation and urinary tract infections*
  - *Ideally suited for long term use*

Führungsdraht, PTFE-beschichtet, flexible Spitze, feste Seele, 150 cm  
 Guide wire, PTFE coated, flexible tip, fixed core, 150 cm




CarboSoft Ureterschiene, beide Enden offen  
 CarboSoft Ureteral stent, both ends open

Schieber, 40 cm  
 Pusher, 40 cm

## CarboSoft Ureteral Stent

DLC coated Soft-Polyurethane

CH	Länge / Length cm		FD-Ø GW-Ø	SET Beide Enden offen Both ends open	Ureterschiene ohne Führungsdraht beide Enden offen Ureteral stent without guide wire both ends open
	A	B	Ø / inch		
6	28	24	.035	3090-1400	3090-1430
	30	26	.035	3090-1500	3090-1530
	32	28	.035	3090-1600	3090-1630
7	28	24	.035	3090-2400	3090-2430
	30	26	.035	3090-2500	3090-2530
	32	28	.035	3090-2600	3090-2630
8	28	24	.038	3090-3400	3090-3430
	30	26	.038	3090-3500	3090-3530
	32	28	.038	3090-3600	3090-3630



Verpackungseinheit: 10 Stück  
 Box: 10 units

Distributed by:

